

Infectious Disease Control III– Level 2

Assignment of values (**REF** BYZKG0210 **LOT** 18C1101)

Instrument	Analyte	Unit	Infectious Disease Control III-Level 2	
			Target	Range
Auto Lumo A2000/Auto Lumo A2000 plus/Auto Lumo A1000/Auto Lumo A6200/AutoLumo S900	HBsAg	IU/mL	0.45	0.27-0.63
	HBeAg	IU/mL	0.32	0.19-0.45
	Anti-TP	S/CO	2.92	1.75-4.08
	Anti-HCV	S/CO	2.65	1.59-3.70
	Anti-HIV (Anti-HIV CLIA Microparticles)	S/CO	2.95	1.77-4.13
	Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles)	S/CO	2.59	1.55-3.63
	Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles) (RU)	S/CO	1.10	0.55-1.65
Auto Lumo A1860	HBsAg	IU/mL	0.44	0.26-0.62
	HBeAg	IU/mL	0.32	0.19-0.45
	Anti-TP	S/CO	2.99	1.79-4.18
	Anti-HCV	S/CO	2.76	1.66-3.87
	Anti-HIV (Anti-HIV CLIA Microparticles)	S/CO	2.72	1.63-3.81
	Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles)	S/CO	2.66	1.60-3.73
	Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles) (RU)	S/CO	1.20	0.60-1.80

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.