

Infectious Disease Control III– Level 4

Assignment of values (**REF** BYZKG0410 **LOT** 11B1001)

Instrument	Analyte	Unit	Infectious Disease Control III-Level 4	
			Target	Range
Auto Lumo A2000/Auto Lumo A2000 plus/Auto Lumo A1000/Auto Lumo A6200/AutoLumo S900	HBsAg	IU/mL	19.23	11.54-26.92
	HBeAg	IU/mL	6.19	3.72-8.67
	Anti-TP	S/CO	11.81	7.08-16.53
	Anti-HCV	S/CO	9.01	5.41-12.61
	Anti-HIV(Anti-HIV CLIA Microparticles)	S/CO	11.96	7.18-16.75
	Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles)	S/CO	11.29	6.77-15.81
Auto Lumo A1860	HBsAg	IU/mL	17.85	10.71-25.00
	HBeAg	IU/mL	6.74	4.04-9.43
	Anti-TP	S/CO	12.69	7.61-17.77
	Anti-HCV	S/CO	10.40	6.24-14.56
	Anti-HIV(Anti-HIV CLIA Microparticles)	S/CO	12.20	7.32-17.08
	Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles)	S/CO	12.45	7.47-17.43

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.