

Infectious Disease Control III– Level 2

Assignment of values (**REF** BYZKG0210 **LOT** 11B1001)

Instrument	Analyte	Unit	Infectious Disease Control III-Level 2	
			Target	Range
Auto Lumo A2000/Auto Lumo A2000 plus/Auto Lumo A1000/Auto Lumo A6200/AutoLumo S900	HBsAg	IU/mL	0.52	0.31-0.73
	HBeAg	IU/mL	0.35	0.21-0.49
	Anti-TP	S/CO	2.89	1.73-4.04
	Anti-HCV	S/CO	2.85	1.71-3.99
	Anti-HIV (Anti-HIV CLIA Microparticles)	S/CO	2.54	1.52-3.55
	Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles)	S/CO	2.49	1.50-3.49
Auto Lumo A1860	HBsAg	IU/mL	0.50	0.30-0.70
	HBeAg	IU/mL	0.37	0.22-0.52
	Anti-TP	S/CO	3.05	1.83-4.27
	Anti-HCV	S/CO	3.03	1.82-4.24
	Anti-HIV (Anti-HIV CLIA Microparticles)	S/CO	2.76	1.66-3.86
	Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles)	S/CO	2.69	1.61-3.77

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.