

## Infectious Disease Control III – Level 1

Assignment of values (REF BYZKG0110 LOT 20E1102)

| Instrument   | Analyte   | Unit  | Infectious Disease Control III-Level 1 |
|--|---|-------|--|
|  |   |       | Range                                  |
| Auto Lumo A2000/Auto Lumo A2000 plus/Auto Lumo A1000/Auto Lumo A6200/AutoLumo S900/Auto Lumo A1860 | HBsAg   | IU/mL | <0.05                                  |
|  | HBeAg   | IU/mL | <0.10                                  |
|  | Anti-TP   | S/CO  | <1.00                                  |
|  | Anti-HCV  | S/CO  | <1.00                                  |
|  | Anti-HIV (Anti-HIV CLIA Microparticles)             | S/CO  | <1.00                                  |
|  | Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles)      | S/CO  | <1.00                                  |
|  | Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles) (RU) | S/CO  | <1.00                                  |

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.