

Liquid Negative Serology Control

Assignment of values (**REF** CFCB0101 **LOT** 11H1101)

Instrument	Analyte	Unit	Liquid Negative Serology Control		
			Target	SD	Range
Auto Lumo A2000/A2000 plus	HBsAg	IU/ml	/	/	<0.05
	HBsAb	IU/ml	/	/	<10
	HBeAg	IU/ml	/	/	<0.1
	HBeAb	IU/ml	/	/	<0.15
	HBcAb	PEIU/ml	/	/	<0.7
	Anti-HIV	S/CO	/	/	<1
	Anti-TP	S/CO	/	/	<1
	Anti-HCV	S/CO	/	/	<1
	Anti-HEV IgM	S/CO	/	/	<1
	Anti-HEV IgG	S/CO	/	/	<1
	Anti-HAV IgM	S/CO	/	/	<1
	Anti-HBc IgM	S/CO	/	/	<1
	HBV-PreS1 Ag	AU/ml	/	/	<6
	HIV Ag/Ab	S/CO	/	/	<1

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.