

Hypertension Markers Control– Level 2

Assignment of values (**REF** CFA0302 **LOT** 30H1101)

Instrument	Analyte	Unit	Hypertension Markers Control Level 2		
			Target	SD	Range
Auto Lumo A2000 plus/A6200/A1860	ACTH	pg/mL	203.58	15.22	157.93-249.23
	AlI	pg/mL	291.18	24.66	217.21-365.14
	ALD	pg/mL	305.64	27.37	223.54-387.75
	Cortisol(IVDD)	µg/dL	26.35	2.16	19.86-32.85
	Cortisol(IVDR)	µg/dL	23.41	1.77	18.11-28.71
	Renin	pg/mL	84.42	7.16	62.95-105.90
	3-MT	pg/mL	84.28	6.96	63.40-105.17
	MN	pg/mL	318.85	25.56	242.17-395.52
	NMN	pg/mL	265.06	20.80	202.68-327.45

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.