

Hypertension Markers Control– Level 1

Assignment of values (**REF** CFA0301 **LOT** 30H1101)

Instrument	Analyte	Unit	Hypertension Markers Control Level 1		
			Target	SD	Range
Auto Lumo A2000 plus/A6200/A1860	ACTH	pg/mL	37.10	3.13	27.70-46.51
	AlI	pg/mL	103.34	7.51	80.82-125.87
	ALD	pg/mL	116.72	9.86	87.12-146.31
	Cortisol(IVDD)	µg/dL	9.38	0.70	7.28-11.48
	Cortisol(IVDR)	µg/dL	8.79	0.68	6.74-10.83
	Renin	pg/mL	15.55	1.47	11.15-19.96
	3-MT	pg/mL	29.58	2.85	21.02-38.14
	MN	pg/mL	91.24	9.52	62.68-119.80
	NMN	pg/mL	78.60	6.99	57.63-99.56

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.